


DQM and Database

Nandita Raha, Antonio Gioiosa

Summary

- Current Status of the DQM
 - ✧ Brief Architecture of the DQM
 - ✧ Source monitor plots – run summary (actual values)
 - ✧ Normalized plots of PMT, PiDs
- Conditional Database
 - ✧ Uses of this database
 - ✧ Future plans

Brief Architecture – Basics of DQM



Overall picture – my contribution

- ✧ Run summary - Source monitor plots (actual values)




- ✧ Normalized plots of PMT, PiDs, and PiD2:PiD1




Run summary plots of SM signals – PiD1, PiD2, PMT


Laser 1 sync amplitude history




Laser 2 sync amplitude history




Laser 3 sync amplitude history




Laser 4 sync amplitude history



Laser 5 sync amplitude history



Laser 6 sync amplitude history





event num

event num

event num


Verification SM signals – Laser 2

Laser 2 sync amplitude history




Overlaid SM signals –PiD2:PiD1, Normalized PiD1, PiD2, PMT


Laser 1 normalized sync amplitude history




Laser 2 normalized sync amplitude history




Laser 3 normalized sync amplitude history




Laser 4 normalized sync amplitude history





Laser 5 normalized sync amplitude history



Laser 6 normalized sync amplitude history



Verification for Laser 2



$$\text{PiD2:PiD1} = 1874.96 / 1687.88$$

1.11084


ConDB – Conditional DataBase

Conditional Database

- ✧ Definition and uses of this database
- ✧ Current work and future plans

ConDB – Conditional DataBase

- Conditions data record the state of the system at the time events were collected.
- Conditions data vary with time. Each value, or set of values, describes the state of the system during a limited interval of time and referred to as the Interval Of Validity (IOV.)



- Useful when many clients run simultaneously on interactive and GRID resources – so read this data is stored in central databases or files.
- Well applicable for data changing with time – like as **calibration**, alignment, attenuation, pedestal, etc. for detector channels.

Laser Calibration from FHICL files to ConDB

Laser OOF calibration data in a zillion FHICL files. Use cat:

```
> cat /cvmfs/gm2.opensciencegrid.org/specials/OOFDB/60h/  
run15974/subrun99.fcl
```

```
> xtalCorr: [ 0.93647, 0.947338, 0.935881, 0.938746, ...1]
```

Array size 1297 (24 calo x 54 crystals + 1 bool)

Transfer these to the conDB using some scripts (by Igor)

Current Status:

Look at data at

[http://dbweb5.fnal.gov:9090/gm2_con_prod/app/data
f=oofgain_condb_test&t=15921000000](http://dbweb5.fnal.gov:9090/gm2_con_prod/app/data
f=oofgain_condb_test&t=15921000000)

t has run and subrun number. 15921000000 => run 15921 subrun 0.

Structure of the data table

The structure of calibration table is currently like:

15921000000.000000

15921000001.000000

channel, corrconst, isgood

bigint, double precision, boolean

100, 0.982577, True

101, 0.986715, True

102, 0.982656, True

....

1st column **channel**: last 2 digits – xtal index (0 to 53); First 1 or 2 digits - calorimeter index.

2nd column **corrconstant**: the gain constant.

3rd column **isgood**: a flag marking if this row is valid. Currently the third column is read from the last number in the 1297 number arrays from the out-of-fill constants FHiCL files stored on cvmfs.

ConDB - Tags

Tagging:

Users can put a text tag on a particular condition, like version “v1.0” on the database.

Users can tag before recalibrating, to remember their current state, and they can tag after recalibration – version control

Supersede a calibration set with a better one and keep it all straight in the DB

Uses and Future Steps

- DQM - helpful for debugging purposes, checking stabilities, comparing with averages etc.
- DQM - Summary histograms for debugging purposes for SM implemented (can add LM?)
- ConDB – only a few subruns have been added - ongoing work in progress
- Plan to tag calibration sets for various run groups in ConDB


Back up

Unstable – SM PiD1 for laser 1, Run 18151??

NUMBER 18151 EVENT 3412 N. LM ISLANDS (for each channel) N. SM ISLANDS (for each channel)


LASER 1

Laser 1 sync amplitude history



LASER 2

Laser 2 sync amplitude history



LASER 3

Laser 3 sync amplitude history

